



## PRECINT SHUTTLE IV Light Curing Unit

- UVA light curing unit for the polymerisation of light curing material and synthetic material
- Big enough to take up to 4 models
- Already preset polymerisation time
- Optical and acoustic signal after the polymerisation process is finished
- 4 light tubes with an emitter zone of 450 nm
- Can also be used for some composites and LC block out material
- One-Touch switching
- Low energy consumption during stand-by operation

**NEW**



Get a FREE Yeti Light Cure Die Block Gel worth £12.95

Polymerisation time	4 min 50 sec	
Power	35 Watt	
Size	Opening size	55mm x 135 mm
	Drawer internal size	130mm x 155mm

Y4949000 **£290**

MESTRA 65th BIRTHDAY

## Plaster Trap

- Made of resistant polypropylene
- Quick and easy bag change
- Air tight seal – prevents smells and spread of bacteria
- Designed to fit under most sink units
- H 35cm, Diameter 35cm
- Capacity 14 litres

**£175** EACH



VISIT OUR WEBSITE WWW.BRACON.CO.UK

## Advanced materials for today's dental laboratory.



### Dip Coping Wax

A dipping wax with the perfect balance of elasticity, memory resilience, and dimensional integrity. Produce perfect dip wax copings with sealed margins. Ash free.

HT5092 Yellow Dip Coping Wax, 284 gr. **£24.95**



### MS1 Master Die Hardener/ Sealer

This high quality die hardener and sealer conditions and protects the gypsum die surface with a zero micron resinous coat. Quick drying and smear free.

HT7175 MS1 Master Hardener/ Sealer, 15 ml **£12.95**  
 HT7176 MS1 Master Hardener/ Sealer, 236 ml **£32.95**  
 HT7177 MS1 Master Hardener/ Sealer, Pen **£4.95**



Available from April 1st

### MS3 Master Die Separator

A magical separator for Wax, Porcelain, and Composite restorations. Evaporates quickly, leaving a non-oily surface with no layer buildup.

HT7182 MS3 Master Die Separator, 15 ml **£12.95**  
 HT7183 MS3 Master Die Separator, 236 ml **£32.95**  
 HT7184 MS3 Master Die Separator, Pen **£4.95**



Now brought to you exclusively by



Made in the U.S.A.

## Wax Sprue System by Harvest.

The Ultimate Solution for Precision Casting of Dental Alloys.



### Reservoir Sprues

Used for direct casting, the scientifically designed 1/4" melting reservoir guarantees optimum alloy flow, and eliminates porosity. 8 gauge.

HT5016 Reservoir, (250) **£12.95**  
 HT5018 Reservoir, (500) **£19.95**



### Feeder Long Sprues

Get a return on your investment with consistent, high quality castings with perfect fit and accurate margins. Nozzle design saves metal. 8 gauge.

HT5011 Feeder Long, (250) **£12.95**  
 HT5012 Feeder Long, (500) **£19.95**



### Sprue Buttons

Facilitate accurate and correct direct sprue architecture. Acts as a funnel to guide molten alloy from the crucible into the pattern void.

HT5019 Buttons, (100) **£8.95**



### Runner Bars

A specially engineered solid center reservoir promotes optimal alloy fluidics resulting in casting completeness. 6 gauge.

HT5015 Runner Bar, (50) **£14.95**  
 HT5017 Runner Bar, (200) **£35.95**



### Circular Sprue

Offering a large number of precision castings at one shot. A special flow-dynamic design allows smooth alloy flow with no turbulence. 6 gauge.

HT5028 Large, 1-9/16 dia (25) **£13.95**  
 HT5029 Large, 1-9/16 dia (100) **£34.95**



### Feeder Short Sprue

Wax sprue gates join patterns to the distribution bar to ensure adequate metal supply. Ensures dense, accurate castings. 8 gauge.

HT5009 Feeder Short (500) **£13.95**

TotalConfidence

www.harvestdental.com

Made in the U.S.A.

# BULLETIN 153



BRACON SOUTHERN: TEL: 01580 817000 FAX: 01580 819455 BRACON NORTHERN: TEL: 0161 480 6731 E-MAIL: SALES@BRACON.CO.UK

## Sarros Sculpturing Waxes

- Our new economic C&B waxes
- Available in *std opaque* and *ash free* types
- Working temperature 65.1°C
- Solidus temperature 55.7°C
- All in 70g tins

FREE SAMPLES NOW AVAILABLE

SCULPTURING		ASH FREE	
SS1013	Grey	SS1018	Grey
SS1014	Beige	SS1019	Beige
SS1015	Blue	SS1020	Blue
SS1016	Green	SS1021	Green
SS1017	White	SS1022	White

All available at only **£12 per tin**



BUY 1 GET 1 FREE

## GENERAL MATERIALS

### Gemini silicone

- *Superb detail* and stable
- Shore hardness of *Gemini Standard* – 8° at 30 min
- *Seriously accurate* low shrinkage of just <0.1%

SP8015 – 2kg pack @ only **£39.95** each

OR BUY:

**2 Packs for £70**



ALL PRICES SUBJECT TO VAT – ALL OFFERS IN THIS LEAFLET END 30 APRIL 2011



### Dr Denti denture boxes

- Durable, leak proof plastic case
- 80ml capacity: enough to soak a F/F
- Includes easy grip cleaning brush and mirror

SG1000

**£4.95** each

BUY  
12 for only **£4.50** each  
48 for only **£3.95** each

### PROSTHETICS MATERIALS & TEETH



### EQUIPMENT

#### Denstar 140

Digitally controlled, single electric wax knife

- Highly accurate temperature control
- Temp range: 50° – 200° C
- Memory function for frequently used temperatures
- Error function – protects the unit from accidental electrical damage
- Magnetic hand-piece holder
- Convenient drawer for hand-piece tips
- Ergonomic, sturdy design



**SAVE OVER 20%**

ED1140 – Only **£110**

### PROSTHETICS MATERIALS & TEETH

#### Acrifix Light curing denture resin gel

- Stable gel – easily and accurately positioned then light cured
- MMA free
- Time saving and sparing in use
- Cures in all standard light boxes

#### Uses:

- Completion of small saddles on partial or implant dentures
- Fabrication of small anterior prosthesis
- Correction of bubbles, pores and cracks on denture base
- Quick repairs of broken off teeth
- Relining of saddles on partial dentures



SS2250 Acrifix pink kit (3 x 3g, bonder and tips) Only **£55**

LEASING OPTIONS ARE NOW AVAILABLE FOR ADVICE PLEASE CALL IAN LEE ON 01580 817001

### ORTHODONTIC MATERIALS



#### Pallarium Solder

- "Cored" (contains flux) solder
- Melting point of 580°C
- Solders nickel chrome, stainless steel, chrome cobalt, gold etc
- Available on white or yellow colours
- Packs of 12 sticks

SS7001 White **£30** each

SS7002 Yellow **£30** each

### GENERAL MATERIALS

#### Megaview Head Loupe

- Comfortable wide field of view and light weight
- Adjustable and can be worn with your own glasses
- 3 lenses different strengths
- Flip-up feature



SK7140 Special offer at just **£85**

VISIT OUR WEBSITE WWW.BRACON.CO.UK

### Eve Diasynt Plus and Diacera

#### Diasynt Plus for removing and shaping

- High performance, synthetically bound diamond abrasives
- Developed specifically for highly resistant materials such as Zirconia
- High removal rate but very low heat generation – minimising the risk of micro-cracks
- Efficient and economical
- Available in two grades: Green coarse and Pink Medium
- Use at 8 – 12,000 rpm

#### Diacera polishers

- Two step polishing system for Zirconia
- Green – Smoothing and pre-polishing
- Orange – High lustre polishing



E19154	DYP13g DIASynt PLUS	EACH	£8.95
E19155	DYP14g DIASynt PLUS	EACH	£9.45
E19156	DYP9g DIASynt PLUS	EACH	£9.45
E19157	DYP15g DIASynt PLUS	EACH	£15.95
E19158	DYP8g DIASynt PLUS	EACH	£9.45
E19164	DYP13m DIASynt PLUS	EACH	£9.45
E19165	DYP14m DIASynt PLUS	EACH	£8.95
E19166	DYP9m DIASynt PLUS	EACH	£9.45
E19167	DYP15m DIASynt PLUS	EACH	£15.95
E19168	DYP8m DIASynt PLUS	EACH	£9.45
E15814	H2DCmf DIACERA	EACH	£7.95
E15816	H8DCmf DIACERA	EACH	£7.95
E15813	H2DC DIACERA	EACH	£7.95
E15815	H8DC DIACERA	EACH	£7.95



### EQUIPMENT



#### Zhermack P03 Hydraulic Press

- Clamps up to 3 flasks
- Equipped with leverage operated pump
- Suitable for flasks up to 124mm wide
- Very robust and durable

EE8305800 **£395**

At Bracon we have 2 full time engineers to service and repair your equipment please call 01580 817000 for more information

VISIT OUR WEBSITE WWW.BRACON.CO.UK

### GENERAL MATERIALS

#### Yeti Flavoured bite waxes

- Stable when warm or cool for excellent bite registration
- Flavoured for patient acceptance
- Colours/flavours correspond to hardness:

Red: Strawberry flavour = Hard  
Green: Apple flavour = Medium  
Yellow: Lemon flavour = Soft

Y6120490	Yellow	490g	<b>£15.95</b>
Y6121850	Yellow	1850g	<b>£59.95</b>
Y6130490	Red	490g	<b>£15.95</b>
Y6131850	Red	1850g	<b>£59.95</b>
Y6140490	Green	490g	<b>£15.95</b>
Y6141850	Green	1850g	<b>£59.95</b>



### EVE Techniks

Three step polishing system for removing, smoothing and polishing of prostheses and orofacial orthopaedic acrylics



#### Removing and shaping

Item:	H334K	H334K	H332K	H332K	H330K	H330K
Item No.:	E13174	E13174S	E13178	E13178S	E13240	E13240S
Quantity:	P100	EACH	P100	EACH	P100	EACH
Price:	£114.95	£1.45	£114.95	£1.45	£114.95	£1.45

#### Smoothing

Item:	H336Km	H336Km	H334Km	H334Km	H332Km	H332Km	H330Km	H330Km
Item No.:	E13157	E13157S	E13175	E13175S	E13177	E13177S	E13250	E13250S
Quantity:	P100	EACH	P100	EACH	P100	EACH	P100	EACH
Price:	£114.95	£1.45	£114.95	£1.45	£114.95	£1.45	£114.95	£1.45

#### High-Lustre Polishing

Item:	H336Kf	H336Kf	H334Kf	H334Kf	H330Kf	H330Kf
Item No.:	E13167	E13167S	E13176	E13176S	E13260	E13260S
Quantity:	P100	EACH	P100	EACH	P100	EACH
Price:	£114.95	£1.45	£114.95	£1.45	£114.95	£1.45

VISIT OUR WEBSITE WWW.BRACON.CO.UK

### Eve Diapol polishers

- 3 – Stage polishing system for ceramic and metal
- Top quality, diamond impregnated polishers
- Very fast, efficient in use
- Extremely high quality and long lasting
- BLUE (coarse) for removal and shaping
- PINK (medium) for smoothing
- GREY (fine) for high lustre polishing



#### Removing and shaping



#### Smoothing



E17226	L26Dg	26 x 2 (mm)	EACH	£26.95
E15201	SL20Dg	20 x 1	EACH	£16.95
E17212	H2Dg	4 x 13	EACH	£7.95
E17208	H8Dg	11 x 2	EACH	£7.95
E17126	L26Dmf	26 x 2	EACH	£26.95
E15002	SL20Dmf	20 x 1	EACH	£16.95
E17112	H2Dmf	4 x 13	EACH	£7.95
E17108	H8Dmf	11 x 2	EACH	£7.95
E15126	L26D	26 x 2	EACH	£26.95
E15203	SL20D	20 x 1	EACH	£16.95
E15102	H2D	4 x 13	EACH	£7.95
E15108	H8D	11 x 2	EACH	£7.95

### GENERAL MATERIALS

#### Bracon Matro tungsten cutters

- Designed specifically for use on acrylics and soft liners
- No "snagging" when new
- Leaves smooth surface for easy finishing
- Range of 6 shapes

315102	<b>£8.60</b>
314902	<b>£8.60</b>
314802	<b>£8.60</b>
312801	<b>£10.44</b>
315002	<b>£8.60</b>
313001	<b>£11.78</b>

Buy any 5 + get the 6th FREE



VISIT OUR WEBSITE WWW.BRACON.CO.UK